

# Home Water Quality Test Kit

Analyte	Procedure	Negative	Positive		
Hardness	Add Pillow RW-1030 and mix	Blue	Red		
Chlorine	Add Pillow RW-230 and mix	Colorless	Purple		
Iron	Add Pillow RW-280 and mix – 5 min	Colorless	Orange		
Nitrates	Add Pillow RW-8310 and mix – 10 min	Colorless	Red		
pH	Dip the pH strip in the water sample and take it out promptly. Find pH by comparing the colors to the chart below				
5.0	6.0	7.0	8.0	9.0	10.0

Have your prospects find what is the water quality problems they want you to correct:

- The kit contains tests for: Hardness, Chlorine (also chloramines react), Iron, Nitrate and pH in an easy to use format.
- You can include your business card and the bag has a hole suitable to hand on door knobs, stair railings, etc.
- Very low cost to create awareness of the water quality in the areas you service



Test Kits Components

(Insert Your Business Card)



Test Kit Hanging on a Door Knob